gampt **ULTRASOUND ● SCHOOL**

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


The screenshot displays the ImageViewSchool software interface. The main window shows a B-mode ultrasound image of a heart. Two measurement lines are drawn on the image: a yellow line labeled '1' and an orange line labeled '2'. A green vertical line is also visible. On the right side, there is a depth scale from 0.0 to 90.0 mm. The software interface includes a menu bar (Exit, Print, About, Reset), a parameter panel with 'Stop' and 'Freeze' buttons, and a control panel with 'TGC', 'A-Line', and 'Presets' tabs. The 'Main' tab is active, showing 'Frequency [MHz]' (2.0, 3.0, 4.0, 5.0) and 'Depth [mm]' (210). The 'Info' panel on the right shows 'Probe: C5-2R60S-3' and 'Cursor' options (Line, Ellipse). The 'Lines' panel shows two lines with lengths: 1 (8.76 mm) and 2 (7.79 mm). The 'State' panel shows 'Imaging is frozen ...'. The logo 'gamPT ULTRASONIC SOLUTIONS 1.0.19.4.2' is visible in the bottom right corner of the software window.

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The screenshot displays the ImageViewSchool software interface. On the left, a hand holds a white ultrasound probe against a pig's ear. The software window shows a B-mode ultrasound image of the ear. A vertical green line is drawn through the image, with two orange lines forming a crosshair. A depth scale on the right side of the image ranges from 0.0 to 90.0 mm. The software interface includes a menu bar (Exit, Print, About, Reset), a Parameter panel with 'Stop' and 'Freeze' buttons, and a main panel with 'TGC', 'A-Line', and 'Presets' tabs. The 'Main' tab is active, showing 'Frequency [MHz]' with options 2.0, 3.0, 4.0, and 5.0 (selected). The 'Info' panel on the right shows 'Probe: C5-2R60S-3' and 'Cursor' options (Line selected, Ellipse unselected). The 'Lines' panel shows two lines with lengths: Line 1 (17,13 mm) and Line 2 (16,96 mm). The 'State' panel shows 'Imaging is frozen ...'. The GAMPT logo is visible in the bottom right corner of the software window.

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The screenshot displays the ImageViewSchool software interface. On the left, a hand holds an ultrasound probe over an orange. The software window shows a B-mode scan of the orange with two yellow lines and orange markers. The parameter panel on the left includes a 'Freeze' button, 'A-Line' and 'Presets' tabs, and a frequency selection of 5.0 MHz. The main scan area shows two yellow lines with orange markers, labeled 1 and 2. The right panel shows a depth scale from 0.0 to 90.0 mm and a green waveform. The info panel on the far right shows 'Probe: C5-2R60S-3', 'Cursor' type (Line selected), and 'Lines' with lengths: 1 (20,86 mm) and 2 (10,25 mm). The state is 'Imaging is frozen ...'. The gamPT logo is visible in the bottom right corner of the software window.

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